

Product Description and Aim of Treatment


Gerbinol® CF is a casein-free powdered adsorbent of tannic substances. Gerbinol® CF smoothens unharmonious wines with components in unbalanced proportions, wines which - since rich in tannin - are rough, lacking softness and an agreeable and clean sensation. Gerbinol® CF reduces, respectively removes also off-flavours. Gerbinol® CF is based on different gelatins, silicates and isinglass.

Permitted according to current laws and regulations of the EU. Purity and quality are proved by specialized laboratories.

Dosage and Application

| | |
|----------------------------------|-----------------|
| Slight sensory correction: | 5 – 10 g/100 L |
| Treatment of tannic substances: | 15 – 30 g/100 L |
| Correction of sensory disorders: | 30 – 50 g/100 L |

Adsorbent of tannic substances, based on gelatins, silicates and isinglass with softening and harmonising effect, casein-free



Principally Gerbinol® CF can be applied in every stage, i.e., also before the usual clarifying treatment. It is however recommendable to perform Gerbinol® CF treatment only after completed clarification, since only then a sensory unbalance or unclear sensory impression is clearly perceptible.

It is recommended to conduct pretests.

After having determined the optimal dosage, Gerbinol® CF is mixed with a tenfold amount of water (if possible, lukewarm) under intensive stirring to assure that the mixture is free of lumps. Afterwards add to the container equally under permanent stirring. Gerbinol® CF is a special product and thus the preparation of the mixture for subsequent wine treatment has to be done properly, with the appropriate care.

Storage

Store in a dry place. Protect from foreign odours and moisture. Reseal opened packagings tightly